



A modular tool to aggregate results from bioinformatics analyses across many samples into a single report.

Report generated on 2022-07-21, 14:00 based on data in: /Users/pp/Documents/H3K27ac/fastq

General Statistics

Copy table

Configure Columns

Plot

Showing 21/21 rows and 4/5 columns.


Sample Name	% Dups	% GC	Read Length	M Seqs
GSM3752845_Y-15T-H3K27ac	31.1%	46%	35 bp	22.7
GSM3752848_Y-13A-H3K27ac	33.6%	47%	28 bp	21.8
GSM3752849_Y-18T-H3K27ac	30.0%	46%	34 bp	25.0
GSM3752850_Y-14A-H3K27ac	33.7%	46%	27 bp	22.2
GSM3752851_Y-19T-H3K27ac	30.8%	45%	34 bp	25.1
GSM3752852_Y-15A-H3K27ac	26.9%	48%	36 bp	23.8
GSM3752854_O-8A-H3K27ac	34.2%	48%	29 bp	22.9
GSM3752855_O-10T-H3K27ac	29.9%	43%	37 bp	22.4
GSM3752856_O-6A-H3K27ac	28.7%	46%	35 bp	23.8
GSM3752857_O-11T-H3K27ac	0.0%	0%	0 bp	0.0
GSM3752858_O-7A-H3K27ac	33.7%	46%	33 bp	24.7
GSM3752861_O-9A-H3K27ac	26.8%	46%	38 bp	21.4
GSM3752862_O-10A-H3K27ac	29.4%	45%	35 bp	24.6
GSM3752863_O-11A-H3K27ac	24.0%	44%	43 bp	20.0
GSM3752864_AD-1T-H3K27ac	35.4%	46%	30 bp	26.3
GSM3752866_AD-2T-H3K27ac	39.5%	47%	28 bp	28.0
GSM3752867_AD-2A-H3K27ac	29.5%	44%	37 bp	22.9
GSM3752868_AD-3T-H3K27ac	26.1%	46%	37 bp	22.6
GSM3752870_AD-3A-H3K27ac	32.5%	44%	34 bp	25.6
GSM3752871_AD-5A-H3K27ac	39.0%	48%	28 bp	28.4
GSM3752872_AD-8T-H3K27ac	33.0%	45%	35 bp	25.0

FastQC

FastQC is a quality control tool for high throughput sequence data, written by Simon Andrews at the Babraham Institute in Cambridge.


Sequence Counts

Sequence counts for each sample. Duplicate read counts are an estimate only.




Sequence Quality Histograms

The mean quality value across each base position in the read.



Per Sequence Quality Scores

The number of reads with average quality scores. Shows if a subset of reads has poor quality.




Per Base Sequence Content

7

The proportion of each base position for which each of the four normal DNA bases has been called.

Click a sample row to see a line plot for that dataset.

Rollover for sample name




Per Sequence GC Content 21

The average GC content of reads. Normal random library typically have a roughly normal distribution of GC content.


Percentages

Counts




Per Base N Content 21

The percentage of base calls at each position for which an N was called.




Sequence Length Distribution 20

The distribution of fragment sizes (read lengths) found. See the FastQC help



Sequence Duplication Levels 10 11


The relative level of duplication found for every sequence.



Overrepresented sequences

11

The total amount of overrepresented sequences found in each library.



Adapter Content

20

The cumulative percentage count of the proportion of your library which has seen each of the adapter sequences at each position.

No samples found with any adapter contamination > 0.1%

Status Checks

Status for each FastQC section showing whether results seem entirely normal (green), slightly abnormal (orange) or very unusual (red).

Sort by highlight

Min: 0

Max: 1

